DATASHEET

THE SINGLE-USE FLEXIBLE DIGITAL CYSTOSCOPE

MAXIMIZE PHYSICIAN EFFICIENCY

Scheduling and workflow flexibility to increase daily patient procedures - including last-minute add-ons

MAXIMIZE STAFF PRODUCTIVITY

Reallocate staff resources from time-consuming reprocessing to patient-focused care

BUDGET OPTIMIZATION


No repairs or reprocessing costs, and limited capital expenditures while increasing procedural capacity


SIMPLIFIED COMPLIANCE

Sterile, patient-ready cystoscope without the complexities of reprocessing and monitoring

PATIENT EXPERIENCE

Expand procedure locations to additional rooms or satellite clinics, and enable image\video review with patients





Developed by the world leader in single-use endoscopy, the aScope 4 Cysto is the newest addition to Ambu's line of endoscopes. When combined with the portable 12.8" HD aView2 Advance monitor, this single-use platform provides an effective solution for flexible cystoscopy procedures across various clinical settings.



To download Instructions For Use for the aScope 4 Cysto, scan the QR code.



SPECIFICATIONS

Optical system	
Field of view	120°
Depth of field	3 - 100 mm
Illumination method	LED
Insertion portion	
Bending angle*	210° up, 120° down
Insertion cord diameter	16.2 Fr / 5.4 mm (0.21")
Distal end diameter	16.2 Fr / 5.4 mm (0.21")
Maximum diameter of insertion portion	18 Fr / 6.0 mm (0.24")
Working length	390 mm (15.4")
Working channel	
Minimum instrument channel width**	6.6 Fr / 2.2 mm (0.086")
Will ill flatti il isti al flerit chariffer wiati	0.0117 2.211111 (0.000)
Operating environment, stora	, ,
	, ,
Operating environment, stora	ge & transportation
Operating environment, stora Temperature - transportation	ge & transportation 10 - 40 °C (50 to 104 °F)
Operating environment, stora Temperature - transportation Temperature - operating environment	ge & transportation 10 - 40 °C (50 to 104 °F) 10 - 40 °C (50 to 104 °F)
Operating environment, stora Temperature - transportation Temperature - operating environment Temperature - storage***	ge & transportation 10 - 40 °C (50 to 104 °F) 10 - 40 °C (50 to 104 °F) 10 - 25 °C (50 to 77 °F)
Operating environment, stora Temperature - transportation Temperature - operating environment Temperature - storage*** Relative humidity	ge & transportation 10 - 40 °C (50 to 104 °F) 10 - 40 °C (50 to 104 °F) 10 - 25 °C (50 to 77 °F) 30 - 85%
Operating environment, stora Temperature - transportation Temperature - operating environment Temperature - storage*** Relative humidity Atmospheric pressure	ge & transportation 10 - 40 °C (50 to 104 °F) 10 - 40 °C (50 to 104 °F) 10 - 25 °C (50 to 77 °F) 30 - 85%
Operating environment, storal Temperature - transportation Temperature - operating environment Temperature - storage*** Relative humidity Atmospheric pressure Sterilization	ge & transportation 10 - 40 °C (50 to 104 °F) 10 - 40 °C (50 to 104 °F) 10 - 25 °C (50 to 77 °F) 30 - 85% 80 - 109 kPa
Operating environment, storal Temperature - transportation Temperature - operating environment Temperature - storage*** Relative humidity Atmospheric pressure Sterilization Method of sterilization	ge & transportation 10 - 40 °C (50 to 104 °F) 10 - 40 °C (50 to 104 °F) 10 - 25 °C (50 to 77 °F) 30 - 85% 80 - 109 kPa
Operating environment, stora Temperature - transportation Temperature - operating environment Temperature - storage*** Relative humidity Atmospheric pressure Sterilization Method of sterilization Material	ge & transportation 10 - 40 °C (50 to 104 °F) 10 - 40 °C (50 to 104 °F) 10 - 25 °C (50 to 77 °F) 30 - 85% 80 - 109 kPa ETO
Operating environment, stora Temperature - transportation Temperature - operating environment Temperature - storage*** Relative humidity Atmospheric pressure Sterilization Method of sterilization Material DEHP phthalates	ge & transportation 10 - 40 °C (50 to 104 °F) 10 - 40 °C (50 to 104 °F) 10 - 25 °C (50 to 77 °F) 30 - 85% 80 - 109 kPa ETO Does not contain any DEHP phthalates



Dimensions	ltem no.	Quantity
Ambu® aScope™ 4 Cysto, Standard Deflection	601001000	5 pcs
Ambu® aScope™ 4 Cysto, Reverse Deflection	600001000	5 pcs
Accessories	Item no.	Quantity
Accessories Stopcock 4-Way. 50 psi (3.5 bar). Swivel Nut, Large Bore****	Item no. K12-11412	Quantity 50 pcs









aScope 4 Cysto offers consistent quality every time



High bending angles of 210°/120°







Ambu

Ambu, Inc. 6230 Old Dobbin Lane Columbia, MD 21045 Tel. 800 262 8462 Fax 800 262 8673 ambuUSA.com

Allows for instillation of fluids and insertion of endoscopic accessories